

SMART



10/04	15/04	20/04	25/04	30/04	35/04	40/04	50/04	10/04	15/04	20/04	25/04	30/04	35/04	40/04	50/04		
04 Taper		Starter		Shaping / Finishing		Finishing				04 Taper		Starter		Shaping / Finishing		Finishing	
Torque: 3.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 3.0 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm	

↑ Excellent Cutting Efficiency

↪ Maximum File Flexibility

EndoArt Smart rotary files have fixed %4 taper allow you to work more safely in narrow canals. The developed special heat treatment technology ensures that the Ni-Ti wire has a much more flexible and durable structure and also turns it into a Gold or Blue color. In addition to its high sharpness efficiency, EndoArt Smart Gold has a higher fracture resistance than standard Ni-Ti files. EndoArt Smart Blue allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and it has a much higher fracture resistance than standard Ni-Ti files. It adapts to the shape of the canal in cases with high slopes.

Smart Kit 15/04-20/04-25/04-30/04-35/04-40/04 19/08-10/04-15/04-20/04-25/04-30/04 19/08-15/04-20/04-25/04-30/04-35/04 19/08-16/02-20/04-25/04-30/04-35/04

Smart Gold Rotary Files

- Rotation Files
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank

Control Memory Heat Treated Wire

Convex Triangle Cross - Section

Smart Blue Rotary Files

- Rotation Files
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank

Wire : Gold
Length : 21mm, 25mm, 31mm
Packaging : 6 Files

Wire : Blue
Length : 21mm, 25mm, 31mm
Packaging : 6 Files



WhatsApp
Get information!