

Endo E class

Features

- Suitable for endodontic procedures
- Push type miniature head
- Auto stop reverse feature
- 15 Memory programs - Mode, Gear ratio, Torque, Speed
- E-Type(SDE-EM24E) Endo Handpiece
- On/off foot switch



Auto Stop Reverse



SDE-EM24E for Endo E Class is the E-Type Motor that can be widely compatible with other contra angles.



Specifications control box only

- Power supply : AC100-120 or 220-240V 50/60 Hz
- Speed : 70 - 1,100 rpm (16:1)
- Torque : 0.5 - 9.9 Ncm (16:1)
- Weight : 910 g
- Dimensions : W175 x D125 x H90(mm)
- Contents : Endo Control Unit, Handpiece, Foot control, Handpiece