

EXPERT



Excellent Cutting Efficiency  
Maximum File's Flexibility



EndoArt Expert files work with a reciprocating motion that allows you to perform root canal treatment with a single file. It includes alternative tapers, %6 and %4 angles for single-file systems, and allows you to prepare at a standard angle along the canal. EndoArt Expert files, which can be used with endomotors capable of reciprocating movement, work with a rotation angle adjustment of 150° left and 30° right. The developed special heat treatment technology makes the Ni-Ti wire a much more flexible and durable structure. EndoArt Expert Gold files have higher sharpness and higher fracture resistance than standard Ni-Ti files. EndoArt Expert Blue allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and compared to standard Ni-Ti files, it has a much higher fracture resistance.

Expert Kit R25-R25-R25-R25-R40-R50 R25-R25-R40-R40-R50-R50 R25-R25-R40-R50 R25-R40-R50

Expert Gold Reciprocation Files

- Reciprocation Files
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank

Control Memory  
Heat Treated Wire

S Section



Expert Blue Reciprocation Files

- Reciprocation Files
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank

Wire : Gold  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

Wire : Blue  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

