

RETREATMENT



EndoArt Retreatment rotary files developed for repeated root canal treatments. Thanks to its high cutting efficiency, it is used in the process of removing the materials in the canals. The system consists of three instruments, D1 - D2 and D3 which are designed for cleaning the canal filling in the (D1) coronal, (D2) middle, and (D3) apical thirds. While D1 has an active tip at the first entry into the root canal, the tips of continuing instruments are not cutting.

Retreatment Kit D1-D1-D2-D2-D3-D3

Retreatment Files

Rotation Files
Specially produced for retreatment cases
Designed to remove filling materials like Gutta Percha
x2 times better Cyclic Fatigue Resistance than conventional Files
Better Cutting Efficiency
Great Fracture Resistance
Lasered Taper and Size information on the shank

Wire : Gold
Length : 16mm, 18mm, 22mm
Packaging : 6 Files

Control Memory
Heat Treated Wire
Convex Triangle
Cross - Section



PATH FILE



Path files are used to create a guide (Glide) path. The use of angled files in combination with hand files before use, reduces the risk of step formation and perforation. The size and taper of the file are 16/02. The ends of the file are rounded to prevent the formation of step and zip formation.

Path Gold Files

Rotation Files
Specially produced for creating glide path
Easily reach apex
Better Flexibility
x2 times better Cyclic Fatigue Resistance than conventional Files
Great Fracture Resistance
Lasered Taper and Size information on the shank

Wire : Gold
Length : 21mm, 25mm, 31mm
Packaging : 6 Files

Control Memory
Heat Treated Wire
Quadrangular



WhatsApp
Get information!