



PEDO



Torque: 2.6 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm	Torque: 2.0 Ncm 300 Rpm



Torque: 1.8 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm	Torque: 1.6 Ncm 300 Rpm

EndoArt Pedo rotary files are produced specially for pediatric cases and has a length of 18 mm. The improved heat treatment technology of EndoArt Pedo ensures that the Ni-Ti wire has a much more flexible and durable structure. EndoArt Pedo files have higher sharpness and fracture resistance than standard Ni-Ti files.

Pedo Kit 15/06-20/04-25/04 15/06-20/04-30/04

Pedo Gold Rotary Files

- Rotation Files
- Specially produced for pedodontic cases
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank

Control Memory
Heat Treated Wire

Convex Triangle
Cross - Section



Pedo Blue Rotary Files

- Rotation Files
- Specially produced for pedodontic cases
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank

Wire : Gold
Length : 18mm
Packaging : 3 Files

Wire : Blue
Length : 18mm
Packaging : 3 Files



WhatsApp
Get information!