



ACTION

SX	S1	S2	F1	F2	F3	F4	F5
Starter	Shaping		Finishing				
Torque: 3.0 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm

Excellent Cutting Efficiency

Maximum Flexibility

SX	S1	S2	F1	F2	F3	F4	F5
Starter	Shaping		Finishing				
Torque: 2.6 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm

EndoArt Action rotary files allow you to move comfortably in the canal thanks to the special structure of its surface, while its variable taper is designed to evacuate debris from the canal. The developed special heat treatment technology ensures that the Ni-Ti wire has a much more flexible and durable structure and also turns it into a Gold or Blue color. EndoArt Action Gold files have a higher sharpness efficiency and higher fracture resistance than standard Ni-Ti files. EndoArt Action Blue files allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and it has a much higher fracture resistance than standard Ni-Ti files. It adapts to the shape of the canal in curvature structures.

Action Kit SX - S1 - S2 - F1 - F2 - F3

Action Gold Rotary Files

- Rotation Files
- Variable Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank

Control Memory
Heat Treated Wire

Convex Triangle
Cross - Section



Action Blue Rotary Files

- Rotation Files
- Variable Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank

Wire : Gold
Length : 21mm, 25mm, 31mm
Packaging : 6 Files

Wire : Blue
Length : 21mm, 25mm, 31mm
Packaging : 6 Files



WhatsApp
Get information!