



RubyComp P Posterior Composite

Rubycomp P is micro hybrid composite system which can be the photocurable applied in the posterior field.

Properties

Non stick to instruments
It has surface hardness is like real enamel
It has high intensity
It has a perfect consistence can be shaped/formed easily
It has high filler rate
It resists to chewing force
There is not discoloration after the polimeration.

Packaging

Refill – 4 gr in a syringe
Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories
Shades – A1 – A2 – A3 – A3,5



RubyFlow Flowable Composite

RubyFlow is flowable light curing radiopaque composite.

Properties

Low shrinkage, low wear
High gloss polishability
Beatiful esthetic restorations
Very low surface roughness
Excellent handling
Non sticky
Very good marginal adaptation

Packaging

Refill – 2 gr in a syringe
Kit 1 - 2 gr in a syringe x 2 pcs.
Kit 2 - 2 gr syringe x 4 pcs.
Kit 3 - 2 gr syringe x 4 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories
Shades – A1 – A2 – A3 – A3,5 – B1 - OA2



ISO 4049

CLASS 2A



WhatsApp
Get information!