



RubyFill Nano Composite

RubyFill Nano is radiopaque nano hybrid composite system which can be the photocurable applied in the anterior and posterior field.

Properties

- Universal Nano Composite
- It has light refraction property like a normal tooth due to the nano fillers in it
- It has high filler rate
- It can be polished high with high degree
- It has high resistance to abrasion
- It provides ease of use

Packaging

- Refill – 4 gr in a syringe
- Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories
- Shades – A1 – A2 – A3 – A3,5 – B1 – B2 – C1 – D3 – OA2 – Extra Bleach

RubyComp Nano Composite

RubyComp Nano is radiopaque nano hybrid composite system which can be the photocurable applied in the anterior and posterior field.

Properties

- Higher Esthetic
- Low shrinkage
- It contains Bis- GMA, UDMA, Bis- EMA resine
- It has high filler rate
- It can be polished with high degree
- It has high resistance to abrasion
- It provides ease of use
- Long durability
- Non sticky

Packaging

- Refill – 4 gr in a syringe
- Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories
- Shades – A1 – A2 – A3 – A3,5 – B1 – B2 – C1 – D3 – OA2 – Extra Bleach



ISO 4049

CLASS 2A



WhatsApp
Get information!