




**Motorized GP Flow**  
Provides consistent flow of GP at desirable rate –  
Precludes voids and makes obturation homogeneous –  
Eliminates hand fatigue that resulted from earlier  
gun type devices

**Intuitive Design**  
Sleek and ergonomic hand-piece –  
Simple push button to dispense GP –  
LED indicator reflects when the device is ready for use –

**Powerful Battery**  
Highly efficient 2600mAh Lithium battery –  
Direct Plug & Play - In case of low battery, -  
use the device directly with adaptor

**Fast-Fill Tips**  
– Autoclavable silver tips for multiple usages  
– Silver tips maintain the thermoplastic temperature of  
GP and offers smooth flow inside canal  
– 360° rotating tip for easy filling at all angles  
– Extra long needle offers convenience and visibility  
– Easy to maintain and replace.

**Comprehensive Display**  
– Three customizable memory programs - T1 to T3  
– Temperature range 100°C - 200°C  
– GP flow rate - Low, Medium, High  
– Reserve GP level indicator



Specifications

Heating Time	20 seconds to 200°C
Charging Time	4 hours to complete charge
Hexagonal Wrench	Use the wrench to bend the tip, rotate the tip and clean the device

