



RubyEtch Etching Gel

RubyEtch is an etching gel containing %37 phosphoric acid. Delivered in 3 ml. and 50 ml syringe with very thin application needles.

Properties

Easily Rinse
Blue coloured
Perfect Surface Retention

Packaging

3 ml syringe
3 ml syringe x 2
3 ml syringe x 4
50 ml syringe



RubyEtch HF Porcelain Etch

RubyEtch HF is 9.5% Buffered Hydrofluoric acid in gel form. Using to improve the mechanical bonding between porcelain and adhesive agent.

Properties

Porcelain Restoration Repair
Buffered Formula
Easily Rinse

Packaging

3 ml syringe
Porcelain Repair Kit - 3 ml syringe + Silane



CLASS 2A



WhatsApp
Get information!