



Ruby Pit and Fissure Sealant

Light curing sealant with flouride releasing for occlusal surfaces. It can be used as preventive products.

Properties

Flouride releasing
Radiopaque
Easy application
Low shirinkage

Packaging

2 gr syringe
2 gr syringe x 2
2 gr syringe x 4



FlexDam Gingiva Barrier

Light curing resin-based gingiva barrier to protect soft tissue.

Properties

Easily removable after treatment
It minimizes heat and soft tissue sensivity during curing
Blue shade
5 second curing

Packaging

2 gr syringe
2 gr syringe x 2
2 gr syringe x 4



CLASS 2A



WhatsApp

Get information!