

RubyPlaton Polishing Discs

Has different kind of discs which has different grins and different sizes. It allows excellent results in polishing. Flexible and easy to use

Properties

14 mm blue, yellow, white, green, brown, pink

- **White** – Extra Fine for Polishing
- **Yellow** – Fine for Finishing
- **Green** – Medium for Countoring
- **Blue** – Coarse for Gross Reduction
- **Brown** – X Coarse for Gross Reduction
- **Pink** – XX Coarse for Gross Reduction

Packaging

80 pcs discs and 1 mandren in a plastic box
240 pcs discs and 4 mandren in a plastic box



RubyPlaton Polishing Strips

Has different kind of strips which has different grin sizes. They are produced combining two different types of grins.

Properties

It allows excellent results in polishing. Flexible and easy to use

- **White** – Extra Fine for Polishing
- **Yellow** – Fine for Finishing
- **Green** – Medium for Countoring
- **Blue** – Coarse for Gross Reduction

Packaging

75 pcs strips in a plastic box
80 pcs discs and 1 mandren + 75 pcs strips in a plastic box
240 pcs discs and 4 mandren + 75 pcs strips in a plastic box

