



RubyLiner with Glass Ionomer

Light curing glass ionomer based cavity liner to protect pulp of the tooth. It is used under direct and indirect restorations, like composites, amalgams, ceramics and another metal restorations.

Properties

Excellent biocompatibility
Releases flouride
Radiopaque
Excellent handling characteristics

Packaging

2 gr syringe
2 gr syringe x 2
2 gr syringe x 4



RubyCal LC Calcium Hydroxide

Light curing calcium hydroxide for pulp capping and liner material under dental filling.

Properties

Direct & indirect pulp capping
Light curing
Single component

Packaging

2 gr syringe
2 gr syringe x 2
2 gr syringe x 4



ISO 4049

CLASS 2A



WhatsApp
Get information!