



EndoArt Sodium Hypochloride

It is %5 sodium hypochloride irrigation liquid for root canal treatment.

Properties

It has surfactant increase penetration to dentinal tubules.
Removes smear layer
Effective canal cleaning
Strong antibacterial effect

Packaging

250 ml liquid



EndoArt Citric Acid

It is %40 citric acid irrigation liquid for root canal treatment.

Properties

Removes smear layer from the dentinal walls
Effective canal cleaning
Cleaning calcium hydroxides from the root canals
Improves the mechanical cleaning of the root canal

Packaging

250 ml liquid



Ruby Iodoform Powder

Iodoform is used for root canal treatment. Can be used for filling of root canals after depulping, chronic periodontitis, wound dressing.

Properties

Pure Iodoform (Pharma Grade)
High PH 12,5

Packaging

30 gr powder



Ruby Gutta Solver

RubyGuttaSolver is used for removing of gutta percha materials during re-treatment of root canals.

Properties

%100 Orange Oil
Non toxic, natural
Easy to wash out

Packaging

10 ml liquid