



RubyZoe Zinc Oxide Cement & Eugenol

Zinc oxide and eugenol based dental cement for temporary filling.

Properties

Temporary Filling
 Root canal sealing
 Root canal filling
 Radiopaque
 Base and lining material under fillings
 Easy to mix

Packaging

45gr powder and 15 ml liquid



RubyZoe Eugenol

Eugenol oil mixing with Zinc oxide powder for temporary filling.

Properties

%100 Eugenol (Pharma Grade)
 Easy to mix

Packaging

15 ml liquid



RubyZoe Zinc Oxide Powder

Zinc oxide powder mixing with eugenol for temporary filling.

Properties

%100 Zinc Oxide (Pharma Grade)
 Base and lining material under fillings
 Temporary Fixing crown and bridges
 Indirect pulp capping
 Easy to mix

Packaging

45gr powder