



Rubydent Curing Light

Rubydent Curing Light offers effective and quality curing performance with ergonomic design and metal head.

- 2 mm resin curing in 1 second
- Long life battery / Over 1300 circles under 10 sec mode
- Stainless steel head
- Replaceable and 360 degree rotating head
- Ergonomic design
- Black or white colour options

Model	: MaxCure7 Power
Input	: 100-240VAC 50Hz/60Hz Power
Output	: 5V DC/1A Wave
length	: 385nm-515nm
Light intensity	: 1000mW/cm ² - 2500mW/cm ²
Battery	: 2600mA

CE Class I



WhatsApp
Get information!